

Surgical Instrument Detergent Enzyme Lubricant Cleaners cut costs. Surgical Instrument Detergent Enzyme Lubricant Cleaners clean fast.

The all-in-ONE Surgical Instrument Detergent Enzyme Lubricant Cleaners deliver the four enzymes necessary for removing all forms of proteinaceous bioburden and the surfactant detergents to clean the surface after the enzymes break down the soil. The all-in-ONE is highly concentrated for fast and cost effective cleaning.

The all-in-ONE Surgical Instrument Detergent Enzyme Lubricant Cleaners are:
non-irritating, 100 % biodegradable, phosphate free, and neutral pH.

The cost cutting all-in-ONE Surgical Instrument Detergent Enzyme Lubricant Cleaners are non-irritating to the nose, eyes, and skin. The all-in-ONE will replace multiple products, speed turnaround times, lower costs, and improve the performance of surgical instruments. The all-in-ONE breaks down all forms of proteinaceous bioburden, cleans the surface, and enhances the passive layer of surgical instruments to prevent corrosion. The all-in-ONE, four Enzyme Detergents (Lipase, Amylase, Carbohydrase, and Protease) break down all forms of bioburden and then rinses clean. Enzyme Detergent Cleaners that do not deliver these 4 enzymes cannot remove all forms of proteinaceous bioburden. The all-in-ONE Detergents deliver an average concentration that is 400% higher than the other enzyme detergents. One gallon of the Surgical Instrument Detergent Enzymes can deliver 512 gallons (1938 liters) of cleaning power. The dilution rate is .25 to .5 ounces per gallon (1.9 ml liters to 3.8 ml per liter) depending on the level of bioburden, encrustation, and stains. The easy FOAM-it Surgical Instrument Detergent Enzyme Lubricant Cleaners lay down a thick blanket of long lasting FOAM to hydrate and remove all forms of proteinaceous debris, hard water deposits, and encrusted stains, for faster cleaning action. When water is added the FOAM dissipates to provide a clear 'free-rinsing' cleaning and conditioning solution.



all-in-ONE Enzyme Detergent Lubricant Cleaners



easy FOAM-it Enzyme Detergent Lubricant Cleaners

The highly concentrated all-in-ONE cleaners are aggressive against your most difficult cleaning challenge and guaranteed to be safe for all types of metal, rigid and flexible scopes, lenses, microsurgical instrumentation and reprocessing equipment.

- significantly reduce the cost of enzyme, detergent, stain removing, and lubricating consumables,
- improve turnaround times for soaking, hand washing, and automated cleaning,
- decrease the inventory level, the number of products used, and the space expended for storing products,
- improve the overall condition, appearance and performance of surgical instruments and scopes,
- combine the detergent and enzyme treatments for automated washers, manual cleaning, and presoaking.

UPS-FS-32 easy FOAM-it Surgical Instrument Detergent 4 Enzyme Lubricant Cleaners, 32 ounce applicators, 12 per case
Ergo UPS1-1 all-in-ONE Surgical Instrument Detergent 4 Enzyme Lubricant Cleaners, 4 gallons with 1 pump per case
Ergo UPS1-5 all-in-ONE Surgical Instrument Detergent 4 Enzyme Lubricant Cleaners, 1 ea. 5 gallon jug
Ergo UPS1-15 all-in-ONE Surgical Instrument Detergent 4 Enzyme Lubricant Cleaners, 1 ea. 15 gallon drum
Ergo UPS1-30 all-in-ONE Surgical Instrument Detergent 4 Enzyme Lubricant Cleaners, 1 ea. 30 gallon drum
Ergo UPS1-55 all-in-ONE Surgical Instrument Detergent 4 Enzyme Lubricant Cleaners, 1 ea. 55 gallon drum

The TOSI Test provides an indication of cleaning performance. The amount of red material, which simulates proteinaceous bioburden, that is cleaned from the stainless steel test strip, indicates the cleaning performance. The Test Strip on the Left shows cleaning results using RUHOF Endozime AW Triple Plus. The Test Strip on the Right shows cleaning results using the all-in-ONE 4 Enzyme Detergent. Your greatest cleaning challenge will be your best proof of the all-in-ONE performance. **Nothing** cleans faster or more effectively than the 4 Enzyme Surgical Instrument Detergents. **Nothing** will lower your cleaning costs as much as the all-in-ONE Surgical Instrument 4 Enzyme Cleaner, **guaranteed**.



[ergo-Logistics](http://www.ergo-Logistics.com)

Customer Service: 1-509-747-5027 CustomerService@ergo-Logistics.com

Admin and Technical Support: 1-509-747-5027 Temple@ergo-Logistics.com

Websites: <http://www.ergo-Logistics.com> www.yourCEBA.com www.surgicalinstrumentcleaner.com