

# all-in-ONE Surgical Instrument Detergent 4 Enzyme Cleaners

The all-in-ONE Enzyme Detergent Cleaner with Lubricant and Conditioners is highly concentrated for fast and cost effective: presoaking, hand washing, and automated cleaning.

The all-in-ONE Surgical Instrument Detergent 4 Enzyme Cleaners are: non-irritating, 100 % biodegradable, phosphate free, and neutral pH.

The cost cutting all-in-ONE Surgical Instrument Detergent 4 Enzyme Cleaners are non-irritating to the nose, eyes, and skin. The all-in-ONE will replace multiple products, speed turnaround times, lower costs, and improve the performance of surgical instruments. The all-in-ONE breaks down all forms of proteinaceous bioburden, cleans the surface, and enhances the passive layer of surgical instruments to prevent corrosion. The all-in-ONE, 4 Enzyme Detergents (Lipase, Amylase, Carbohydrase, and Protease) break down all forms of bioburden and then rinses clean. Enzyme Detergent Cleaners that do not deliver these 4 enzymes cannot remove all forms of proteinaceous bioburden. The all-in-ONE Detergents deliver an average concentration that is 400% higher than the other enzyme detergents. One gallon of the Surgical Instrument Detergent Enzymes can deliver 512 gallons (1938 liters) of cleaning power. The dilution rate is .25 to .5 ounces per gallon (1.9 ml liters to 3.8 ml per liter) depending on the level of bioburden, encrustation, and stains. The easy FOAM-it, all-in-ONE Surgical Instrument Detergent 4 Enzyme Cleaner, lays down a thick blanket of long lasting FOAM to hydrate and remove all forms of proteinaceous debris, hard water deposits, and encrusted stains, for faster cleaning action. When water is added the FOAM dissipates to provide a clear 'free-rinsing' cleaning and conditioning solution. The all-in-ONE Enzyme Detergent comes pre-loaded into Instrument and Endoscope cleaning Sponges offering an additional Margin of Safety and high quality cleaning.



all-in-ONE Enzyme Detergent Cleaners



all-in-ONE easy FOAM-it,

The highly concentrated all-in-ONE cleaners are aggressive against your most difficult cleaning challenge and guaranteed to be safe for all types of metal, rigid and flexible scopes, lenses, microsurgical instrumentation and reprocessing equipment.

- significantly reduce the cost of enzyme, detergent, stain removing, and lubricating consumables,
- improve turnaround times for soaking, hand washing, and automated cleaning,
- decrease the inventory level, the number of products used, and the space expended for storing products,
- improve the overall condition, appearance and performance of surgical instruments and scopes,
- combine the detergent and enzyme treatments for automated washers, manual cleaning, and presoaking.

UPS-FS-32 easy FOAM-it, 32 ounce foam applicators: 12 per case

Ergo UPS1-1 all-in-ONE Surgical Instrument Detergent 4 Enzyme Cleaners, 4 gallons with 1 pump per case

Ergo UPS1-5 all-in-ONE Surgical Instrument Detergent 4 Enzyme Cleaners, 1 ea. 5 gallon jug

Ergo UPS1-15 all-in-ONE Surgical Instrument Detergent 4 Enzyme Cleaners, 1 ea. 15 gallon drum

Ergo UPS1-30 all-in-ONE Surgical Instrument Detergent 4 Enzyme Cleaners, 1 ea. 30 gallon drum

Ergo UPS1-55 all-in-ONE Surgical Instrument Detergent 4 Enzyme Cleaners, 1 ea. 55 gallon drum

The TOSI Test provides an indication of cleaning performance. The amount of red material, which simulates proteinaceous bioburden, that is cleaned from the stainless steel test strip, indicates the cleaning performance. The Test Strip on the Left shows cleaning results using RUHOF Endozime AW Triple Plus. The Test Strip on the Right shows cleaning results using the all-in-ONE 4 Enzyme Detergent. Your greatest cleaning challenge will be your best proof of the all-in-ONE performance. **Nothing** cleans faster or more effectively than the 4 Enzyme Surgical Instrument Detergents. **Nothing** will lower your cleaning costs as much as the all-in-ONE Surgical Instrument 4 Enzyme Cleaner, and **we guarantee it.**



[ergo-Logistics](#)

Customer Service: 1-509-747-5027, email: [Support@ergo-Logistics.com](mailto:Support@ergo-Logistics.com)

Technical Support: 1-509-981-1553, email: [Temple@ergo-Logistics.com](mailto:Temple@ergo-Logistics.com)

Websites are available for additional information: [www.yourCEBA.com](http://www.yourCEBA.com) [www.surgicalinstrumentcleaner.com](http://www.surgicalinstrumentcleaner.com)