

yourCEBA

Customer Service: Beth Crafton RN

509-994-2801

Fax: 509-443-1781

[email: Beth@yourCEBA.com](mailto:Beth@yourCEBA.com)

Technical Support: John Temple

1-509-747-5027

[email: Info@yourCEBA.com](mailto:Info@yourCEBA.com)

www.yourCEBA.com

www.surgical-instrument-washer.com



DS50 STEELCO Surgical Instrument Washer Series

External: 25" x 22" x 21" [632x562x535mm]

Internal: 15.75" x 18.9" x 9.65" [400x480x245mm]

The Instrument Washer Disinfector Series DS 50 is ideal for doctor offices and small surgery centers to wash, disinfect and dry surgical instruments.

DS50 Washers are Easy to install, requires no dedicated drain. counter or bench, as single door or pass through.

Thermal disinfection with a temperature of up to 200°F [93°]

Integrated hot air dryer available as option.

EN ISO 15883 HTM 2030 Compliant for washer disinfectors, CE Classified under medical device and FDA registered, UL - listed

DS50 Surgical Instrument Washer

Standard Key features

RS232 interface for printer- or PC connection

LED/LCD operator panel with 3 adjustable cleaning cycles

Build in water softener

Washing and thermal disinfecting temperature adjustable up to 200°F [93°C]

Temperature monitoring by two independent PT-1000 probes

Exhaust steam condenser is standard

Heavy duty wash pump with high flow rate and low pressure for efficient and gentle cleaning.

Chamber made of stainless steel 316L

High efficient triple water filtering

One or two liquid dosing pump(s) standard included

DS50 Surgical Instrument Washer r DS50

Dimensions: (W x D x H)

External: 19.7" x 22" x 21" [500x560x535mm]

Internal: 15.75" x 18.9" x 9.65" [400x480x245mm]

Max. basket capacity: 12.5gal. [47liter]

Electrical connection: 1x115V/60Hz

Single Door

DS50 DRS Instrument Washer Disinfector DS50DRS

Dimensions: (W x D x H)

External: 22" x 22.63" x 21" [560x575x535mm]

Internal: 15.75" x 18.9" x 9.65" [400x480x245mm]

Max. basket capacity: 12.5gal. [47liter]

Electrical connection: 1x115 or 220V/60Hz

Single door

Integrated dry

DS50/2 DRS Surgical Instrument Washer DS50/2 DRS

Dimensions: (W x D x H)

External: 25" x 22" x 21" [632x562x535mm]

Internal: 15.75" x 18.9" x 9.65" [400x480x245mm]

Max. basket capacity: 12.5gal. [47liter]

Electrical connection: 1x115 or 220V/60Hz

Double door pass through

Integrated dryer



yourCEBA

Customer Service: Beth Crafton RN

509-994-2801

Fax: 509-443-1781

email: Beth@yourCEBA.com

Technical Support: John Temple

1-509-747-5027

email: Info@yourCEBA.com

www.yourCEBA.com

www.surgical-instrument-washer.com



Undercounter Surgical Instrument Washers, Countertop Surgical Instrument Washers,

DS500 Washer Racks. with luer lock fittings and silicone adapter connections, clean cannulated surgical instruments **inside and out** eliminating the costs and risks associated with hand washing surgical instruments

DS500 STEELCO Surgical Instrument Washer

External: 23.6" x 24.8" x 33.1" [600x630x840mm]

Internal: 20.5" x 19.7" x 20.9" [520x500x530mm]

1 Manual door

Easy to install, requires no dedicated drain.

Thermal disinfection with a temperature of up to 200°F [93°]

Integrated hot air dryer available

EN ISO 15883 HTM 2030 Compliant to standards for washer disinfectors

Classified under CE medical device and FDA registered

UL – listed

DS500 Surgical Instrument Washer

Standard Key features

Up to 8 instrument trays 20" x 10" x 2" [500x250x50mm]

RS232 interface for printer- or PC connection

LCD operator panel with 40 adjustable cleaning cycles

Build in water softener

Glass viewing window available

Washing and thermal disinfecting temperature adjustable up to 200°F [93°C]

Temperature monitoring by two independent PT-1000 probes

Exhaust steam condenser standard

Heavy duty wash pump with high flow rate and low pressure for efficient and gentle cleaning

Chamber made of stainless steel 316L

High efficient triple water filtering

Two liquid dosing pumps standard included (additional 2 available as option)

decoTwin® cleaning process ready

Surgical Instrument Washer Disinfector DS500C

DS500 C Dimensions: (W x D x H)

DS500 C External: 23.6" x 24.8" x 33.1" [600x630x840mm]

DS500 C Internal: 20.5" x 19.7" x 20.9" [520x500x530mm]

Max. basket capacity: 37gal. [140liter]

Electrical connection: 1x220V, 3x208V/60Hz

**Single Door Undercounter Countertop Surgical Instrument Washers,
Integrated hot air dryer**

Surgical Instrument Washer Disinfector DS500 DRS

DS500 DRS Dimensions: (W x D x H)

DS500 DRS External: 35.5" x 24.8" x 33.1" [900x630x840mm]

DS500 DRS Internal: 20.5" x 19.7" x 20.9" [520x500x530mm]

Max. basket capacity: 37gal. [140liter]

Electrical connection: 1x220V, 3x208V/60Hz

**Single Door Undercounter Countertop Surgical Instrument Washers,
Integrated hot air dryer**

Integrated detergent storage compartment (adds 11.9 inches to width

